

## Sample Requirements for "*Run*"

### 1. Sequencing Reaction

Our system is set up to cope with sequencing reactions of varying quality and intensity. However, the best cost/quality ratio can be obtained by setting up reactions in a total volume of 10  $\mu$ l with 0.5  $\mu$ l of BigDye Terminator Sequencing Mix (either v1.1 or v3.1). The "*Run*" service is available only for full 96-well plates.

### 2. Reaction purification

For best results and to prolong capillary lifetime small molecules (especially salts) and unincorporated dye terminators have to be removed completely prior to electrophoresis. We recommend an SDS/heat treatment to minimize dye blobs followed by purification via Sephadex columns.

### 3. Job submission

Sign up for your job using the "*Run*" service with either BigDye v1.1 or v3.1 on the sequencing homepage (<http://www.genetik.biologie.uni-muenchen.de/sequencing>). The 96 samples can either be named individually or sample names can be uploaded as a tab-delimited text file for the whole plate at once. For this, arrange sample names in Excel in a 12x8 grid and save this grid accordingly.

### 4. Sample delivery

Deposit your **clearly labeled** samples in the fridge in room G03.031 (Sequencing Service im LMU Biozentrum, Großhaderner Str. 2-4, 82152 Martinsried). If deposited by 10:00 AM they will be processed on the same day. The samples have to be in a half-skirted 96-well plate (can be obtained at the Genomics Service Unit) in a total volume of at least 20  $\mu$ l of which at least 10  $\mu$ l should be formamide to stabilize the reaction products. Usually, sequence data can be downloaded from the sequencing homepage on the same day.

### 5. Troubleshooting

The majority of sequencing problems are due to either incorrect template or primer concentration or contaminants in the template. Contaminants known to interfere with the sequencing reaction are: salts (NaCl, NaAc, KAc, KCl), chelators (EDTA, EGTA), proteins, detergents (SDS, Triton X-100), RNA, chromosomal DNA, organic chemicals (ethanol, chloroform, phenol), divalent cations ( $Mg^{2+}$ ,  $Ca^{2+}$ ,  $Mn^{2+}$ ), and excess PCR primers, dNTPs, enzyme, and buffer components from PCR. Be sure to clean your DNA template and check the quantity and quality before sequencing.